



Leica EM RAPID

Pharmaceutical Milling System for Pills


Leica EM RAPID for Pharmaceutical

The Leica EM RAPID milling system has been designed for preparation of pills and tablets prior to investigation with NIR spectrophotometry for Research on Active Pharmaceutical Ingredient Dispersal. With a tungsten carbide or diamond milling tool, pills can be decapsulated and milled without any smearing effect caused by the pills' coating.


The operating concept of the Leica EM RAPID offers versatile features which are directly converted into user benefits to achieve maximum ease of use and optimal results.



Perpendicular view showing graticule for distance definition – optional



Processing observation



Front face observation

Leica Design by W. Hölbl

al Surface Preparation

- Decapsulation of pills without smearing effect
- Adjustable cutting speed between 300–20.000 rpm
- Defined layer removal in steps of 0.5, 1, 10 or 100 μm
- Brilliant ring LED top light and optimised back light illumination
- Low noise extraction and filtration unit with Hepa-filter




Image a




Image b




Image c

Visible (a), derived single component IR image (b) from tablet surface and (c) composite image from 3 major components
Courtesy of Perkin Elmer

The Leica EM RAPID is Fast

A specially designed pill holder allows a fast transfer to the NIR spectrophotometer after milling without removing the sample from the carrier. It can then be reinserted into the Leica EM RAPID for further defined layer removal, thus obtaining accurate and fast layer information. This way no information is lost and the time to the NIR microscope is very fast.

Single pill holder

The pill is mounted on the insert and clamped in the adaptor. After milling the insert can be placed directly in the NIR microscope. The insert dimensions are compatible with the Perkin Elmer tablet Autosampler.



Multiple pill holder

For preparation and investigation of several tablets at once, a smart designed glass slide holder increases throughput considerably, allowing more than one pill to be milled simultaneously.



... is Safe

The milling tool and the sample are located under a protective cover. The instrument will not run until the cover is closed or stops quickly if the cover is opened during milling. The miller ON/OFF switch is integrated into the control lever for the milling action. Hazardous dust of API material is drawn into a low-noise extraction unit fitted with a Hepa-filter. This filtration system is silent and safe.



Control lever with integrated ON/OFF switch

Exhaust tube connection to the extraction and filtration unit



Low-noise extraction and filtration unit fitted with a Hepa-filter.

... is Easy to Use

Integrated viewing system

With the specimen pivot arm the sample can be observed during preparation at an angle between 0° and 60° , directly onto the front face and 90° to the front face for distance determination with the eyepiece graticule (optional). The Leica EM RAPID features brilliant LED top light ring-illumination and optimized backlighting for excellent viewing. By moving the control lever up and down the milling action is performed with the pre-selected step advance of the miller using the ergonomic positioned hand wheel.

Leica EM RAPID ergonomic MZ stereo outfit


The Leica EM RAPID can be ergonomically configured with an MZ stereo microscope in conjunction with the ergo-wedge. The use of the MZ stereo also allows the use of additional equipment such as TV camera and distance definition equipment.



... is Accurate

The Leica EM RAPID is operated by a built-in control panel which provides all necessary functions to achieve highly accurate results.

- ① Milling speed (rpm) on LCD display
- ② Miller total feed advance on LCD display
- ③ Variable cutting speed between 300-20000 rpm allows settings to be adjusted for fragile pills.
- ④ Clear total feed advance to set starting point for defined layer removal of the pill
- ⑤ Hold advance indication while the miller has been retracted without losing milling position information
- ⑥ ON/OFF ring LED illuminator
- ⑦ The Auto-Feed function allows safe and fast layer removal without operating the feed hand wheel. Two switch positions are selectable for single or multiple pill holder.
- ⑧ Fast forward movement of the miller for quick miller approach towards the sample
- ⑨ Fast retraction of the miller and automatic movement to the home position for miller exchange
- ⑩ Pre-selection of the miller advance in steps of 0.5, 1, 10 or 100 μm at each click stop of the hand wheel encoder. Investigations on a defined layer can be accurately performed by stepwise removal of the tablet material.



Milling tools

A selection of design-optimized tungsten carbide and diamond milling tools with a diameter of 8 mm and 12 mm provide high surface quality without any smearing effects caused by the pills' coating.

Leica Microsystems – the brand for outstanding products

Leica Microsystems' mission is to be the world's first-choice provider of innovative solutions to our customers' needs for vision, measurement and analysis of microstructures.

Leica, the leading brand for microscopes and scientific instruments, developed from five brand names, all with a long tradition: Wild, Leitz, Reichert, Jung and Cambridge Instruments. Yet Leica symbolizes innovation as well as tradition.

Leica Microsystems – an international company with a strong network of customer services

Australia:	North Ryde	Tel. +61 2 8870 3500	Fax +61 2 9878 1055
Austria:	Vienna	Tel. +43 1 486 80 50 0	Fax +43 1 486 80 50 30
Canada:	Richmond Hill/Ontario	Tel. +1 905 762 2000	Fax +1 905 762 8937
Denmark:	Herlev	Tel. +45 4454 0101	Fax +45 4454 0111
France:	Rueil-Malmaison	Tel. +33 1 47 32 85 85	Fax +33 1 47 32 85 86
Germany:	Bensheim	Tel. +49 6251 136 0	Fax +49 6251 136 155
Italy:	Milan	Tel. +39 0257 486.1	Fax +39 0257 40 3475
Japan:	Tokyo	Tel. + 81 3 5421 2800	Fax +81 3 5421 2896
Korea:	Seoul	Tel. +82 2 514 65 43	Fax +82 2 514 65 48
Netherlands:	Rijswijk	Tel. +31 70 4132 100	Fax +31 70 4132 109
People's Rep. of China:	Hong Kong	Tel. +852 2564 6699	Fax +852 2564 4163
Portugal:	Lisbon	Tel. +351 21 388 9112	Fax +351 21 385 4668
Singapore		Tel. +65 6779 7823	Fax +65 6773 0628
Spain:	Barcelona	Tel. +34 93 494 95 30	Fax +34 93 494 95 32
Sweden:	Kista	Tel. +46 8 625 45 45	Fax +46 8 625 45 10
Switzerland:	Heerbrugg	Tel. +41 71 726 34 34	Fax +41 71 726 34 44
United Kingdom:	Milton Keynes	Tel. +44 1908 246 246	Fax +44 1908 609 992
USA:	Bannockburn/Illinois	Tel. +1 847 405 0123	Fax +1 847 405 0164

and representatives of Leica Microsystems
in more than 100 countries.

The companies of the Leica Microsystems Group operate internationally in three business segments, where we rank with the market leaders.

● Microscopy Systems

Our expertise in microscopy is the basis for all our solutions for visualization, measurement and analysis of microstructures in life sciences and industry. With confocal laser technology and image analysis systems, we provide three-dimensional viewing facilities and offer new solutions for cytogenetics, pathology and materials sciences.

● Specimen Preparation

We provide comprehensive systems and services for clinical histo- and cytopathology applications, biomedical research and industrial quality assurance. Our product range includes instruments, systems and consumables for tissue infiltration and embedding, microtomes and cryostats as well as automated stainers and coverslippers.

● Medical Equipment

Innovative technologies in our surgical microscopes offer new therapeutic approaches in microsurgery.