



# Leica EM TP

Tissue Processor for Electron and Light Microscopy  
Resin Processing

Living up to Life

*Leica*  
MICROSYSTEMS

# Introducing the Leica EM TP Resin Processor for EM and LM

With 125 years experience in the manufacture of scientific instruments, Leica has developed a new resin tissue processor, the Leica EM TP. From the wealth of knowledge gained in the fields of histology and electron microscopy the EM TP meets all the needs of today's demanding laboratory. This versatile, very compact instrument can be used for resin processing for electron microscopy (EM) and light microscopy (LM). Reproducibility, the need for a low hazard environment and ease of use were the key reasons for developing the EM TP.

The Leica EM TP consists of a base unit onto which an EM or LM processing outfit can be added. Interchange between EM and LM processing takes only moments.

The latest technology is used for interfacing user and instrument. A simple to use control panel allows data storage of up to 99 programmes, all of which can be named. Reagents can be selected from a reagent list which can be customised by the user. If a mistake is made during programming or selection of a field then the Leica EM TP gives an audible warning and shows on the screen what the error means. Processing can be started immediately or with delay. Programming can also be carried out while the instrument is processing.



Designed by Werner Hölbl

Independent vial seals for each vial provide a stable environment for each reagent. The vials sit in a carousel designed for both EM and LM vials. To provide a safe environment for the user in case of reagent spillage the processing chamber is enclosed. The chamber is also ventilated via a fume extraction system.

The specimens are supported in various baskets depending upon their size, the mesh varying in accordance with sample size. The baskets are located on an arm which can move vertically up and down to allow agitation when required. The speed of agitation is variable.

For processing at specific temperatures the Leica EM TP has a heater and cooler system allowing temperature control between +4°C and +60°C even during a delayed start or finish. For EM the system also has pre-heat and pre-cool. This means that all reagents coming into contact with the tissue will be at the desired temperature, important for high quality EM processing.


Programming is via a membrane-covered keyboard. A delay and 'time to end' feature makes setting delay times simple.

The programmes are protected by battery back-up and for specimen safety the instrument has a second battery back-up enabling processing to continue through a power shortage. The heater/cooler (H/C) unit closes down but the carousel continues to operate thus preventing the specimens remaining in processing reagents for too long a period and damaging them.

### Display

(The example shows the screen during programming an EM run)

- Programme number and name
- Vial number and reagent name
- Processing parameter for selected vial
- Total processing time for EM samples




PROG01 EPON  
VIAL01 Glutaraldehyd  
AG1 01h01 20°C CONT  
EM TOTAL 4h06



Selection and adjusting of programmes

Keys for starting and interrupting a programme



*Leica*

Display a programme  
Lock function  
Time setting to adjust clock  
Printing – to print programmes  
and reagent list

Keys for manual operation

- Rotation of carousel
- Lowering carousel
- Rising carousel
- Agitation
- Temperature setting




## EM Processing

24 polypropylene vials of 20 ml are available to the user for one processing run. The vials are attached to the carousel by a twist fit only.

Reagents can be pre-loaded into the vials attached to the carousel in a fume cupboard. Toxic reagents can be sealed with individual vial caps before taking them over to the Leica EM TP. The heater/cooler (H/C) unit flips back allowing direct access of the carousel to the processing area. The H/C unit has pre-heat and pre-cool thus allowing reagents to be at the correct temperature before coming into contact with the tissue.



H/C unit shown in loading position for safe and convenient access



Processing of EM samples



Loading plate and consumable

## LM Processing

12, 100 ml vials can be twist fitted to the carousel for processing large samples or high volume of specimens into resin. The Leica EM TP has a holder for standard histology cassettes thus allowing very large specimens such as bone to be decalcified and processed. CellSafe holders can also be used for processing.



Moveable H/C unit for safe and easy access of carousel



Processing of LM samples



Cassette holder

