



-150°C Cryo Freezer DW-150W200



Key Features

- Alternative for LN2 Storage
- Convenience
- Reliable Performance

More Detail:

http://www.haiermedical.com/products/preservation/150_cryo_freezer/dw-150w200.html

Specifications

Technical Data	
Type	-150 Cryo Freeer
Cabinet Type	Chest
Climate Class	N
Cooling Type	Forced air cooling
Defrost Mode	Manual
Refrigerant	CFC-Free
Performance	
Cooling performance()	-150
Temperature Range()	-126~-150
Control	
Controller	Microprocessor
Display	LED
Electrical Data	
Power Supply(V/Hz)	380/50
Power(W)	2800

Electrical Current(A)	9
Dimensions	
Capacity(L / Cu.Ft)	200/7.1
Net /Gross Weight(approx)	315/385(kg) 694.5/848.8(lbs)
Interior Dimensions(W*D*H)	667*462*650(mm) 26.3*18.2*25.6(in)
Exterior Dimensions(W*D*H)	1570*850*1060(mm) 61.8*33.5*41.7(in)
Packing Dimensions(W*D*H)	1710*985*1195(mm) 67.3*38.8*47.0(in)
Container load (20' /40' /40'H)	6/13/13
Functions	
High /Low Temperature	Y
Remote Alarm	Y
Hot Condenser	Y
Power Failure	Y
High /Low Voltage	Y
Sensor Error	Y
Low Battery	Y
High Ambient Temp	Y
Door Ajar	N/A
Accessories	
Caster	Y
Foot	Y
Test Hole	Y
Shelves / Inner Doors	-/2
Temperature Recorder	Optional
USB Interface	N/A
Rs232 / 485 Interface	Optional
CO2 Backup System	Optional
LN2 Backup System	Optional
Racks And Boxes	Optional

Racks & Boxes(Optional)

Rack Model	Rack Description	Rack Dimensions (W*D*H)	Q'ty of Racks	Q'ty of Boxes	Q'ty of Samples 9*9 box / 10*10 box
DCJ-11-A	2 inch stainless steel rack 11*1	138*142*633(mm) 5.4*5.6*24.9(in)	12	132	10692/13200

Others

Certificate	/
-------------	---