



### Key Features

- Drawer Structure
- Strong Refrigeration
- Two Alarm Functions

### More Detail:

[http://www.haiermedical.com/products/preservation/40\\_deep\\_freezer/dw-40l188.html](http://www.haiermedical.com/products/preservation/40_deep_freezer/dw-40l188.html)

## Specifications

### Technical Data

Type	-40 Deep Freezer
Cabinet Type	Upright
Climate Class	N
Cooling Type	Direct cooling
Defrost Mode	Manual
Refrigerant	CFC-Free
Noise Level (dB)	43

### Performance

Cooling performance( )	-40
Temperature Range( )	-10~-40

### Control

Controller	Microprocessor
Display	LED

### Electrical Data

Power Supply(V/Hz)	220-240/50
--------------------	------------

<b>Power(W)</b>	290
<b>Electrical Current(A)</b>	2.5
<b>Dimensions</b>	
<b>Capacity(L/Cu.Ft)</b>	188/6.6
<b>Net/Gross Weight(approx)</b>	63/71(kg) 138.9/156.5(lbs)
<b>Interior Dimensions(W*D*H)</b>	380*420*1225(mm) 15.0*16.5*48.2(in)
<b>Exterior Dimensions(W*D*H)</b>	597*610*1475(mm) 23.5*24.0*58.1(in)
<b>Packing Dimensions(W*D*H)</b>	659*699*1548(mm) 25.9*27.5*60.9(in)
<b>Container load (20'/40'/40'H)</b>	27/54/54
<b>Functions</b>	
<b>High/Low Temperature</b>	Y
<b>Sensor Error</b>	Y
<b>Accessories</b>	
<b>Caster</b>	Y
<b>Foot</b>	Y
<b>Test Hole</b>	N/A
<b>Drawers</b>	7
<b>Temperature Recorder</b>	N/A
<b>USB Interface</b>	N/A
<b>Others</b>	
<b>Certificate</b>	CE