



Key Features

- Drawer Structure
- Strong Refrigeration
- Two Alarm Functions

More Detail:

http://www.haiermedical.com/products/preservation/40_deep_freezer/dw-40l262.html

Specifications

Technical Data

| | |
|------------------|-----------------|
| Type | -40 Deep Freeer |
| Cabinet Type | Chest |
| Climate Class | N |
| Cooling Type | Direct cooling |
| Defrost Mode | Manual |
| Refrigerant | CFC-Free |
| Noise Lever (dB) | 44 |

Performance

| | |
|------------------------|---------|
| Cooling performance() | -40 |
| Temperature Range() | -10~-40 |

Control

| | |
|------------|----------------|
| Controller | Microprocessor |
| Display | LED |

Electrical Data

| | |
|--------------------|-----------------------|
| Power Supply(V/Hz) | 220-240/50/208-230/60 |
|--------------------|-----------------------|

| | | |
|--|--|-----|
| Power(W) | 340 | 340 |
| Electrical Current(A) | 3.2 | 3.2 |
| Dimensions | | |
| Capacity(L/Cu.Ft) | 262/9.3 | |
| Net/Gross Weight(approx) | 88/93(kg) 194.0/205.0(lbs) | |
| Interior Dimensions(W*D*H) | 480*465*1430(mm) 18.9*18.3*56.3(in) | |
| Exterior Dimensions(W*D*H) | 657*685*1665(mm) 25.9*27.0*65.6(in) | |
| Packing Dimensions(W*D*H) | 732*745*1756(mm) 28.8*29.3*69.1(in) | |
| Container load (20'/40'/40'H) | 24/48/48 | |
| Functions | | |
| High/Low Temperature | Y | |
| Sensor Error | Y | |
| Accessories | | |
| Caster | Y | |
| Foot | Y | |
| Test Hole | N/A | |
| Drawers | 7 | |
| Temperature Recorder | N/A | |
| USB Interface | N/A | |
| Others | | |
| Certificate | CE | |