



Key Features

- Strong Refrigeration
- Five Alarm Functions
- Drawer Structure

More Detail:

http://www.haiermedical.com/products/preservation/40_deep_freezer/dw-40l348.html

Specifications

Technical Data

Type	-40 Deep Freezer
Cabinet Type	Upright
Climate Class	N
Cooling Type	Direct cooling
Defrost Mode	Manual
Refrigerant	CFC-Free
Noise Level (dB)	45

Performance

Cooling performance()	-40
Temperature Range()	-20~-40

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220-240/50
--------------------	------------

Power(W)	360
Electrical Current(A)	2.9
Dimensions	
Capacity(L/Cu.Ft)	348/12.3
Net/Gross Weight(approx)	137/145(kg) 302.0/319.7(lbs)
Interior Dimensions(W*D*H)	535*607*1228(mm) 21.5*23.8*48.2(in)
Exterior Dimensions(W*D*H)	750*840*1860(mm) 29.3*33.1*73.2(in)
Packing Dimensions(W*D*H)	815*945*20408(mm) 32.1*37.2*80.3(in)
Container load (20'/40'/40'H)	14/28/28
Functions	
High/Low Temperature	Y
Remote Alarm	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Accessories	
Caster	Y
Foot	Y
Test Hole	Y
Drawers	12
Temperature Recorder	Optional
USB Interface	Y
Others	
Certificate	CE