



### Key Features

- Five Alarm Functions
- Drawer Structure
- Strong Refrigeration

### More Detail:

[http://www.haiermedical.com/products/preservation/40\\_deep\\_freezer/dw-40l508.html](http://www.haiermedical.com/products/preservation/40_deep_freezer/dw-40l508.html)

## Specifications

### Technical Data

|                  |                  |
|------------------|------------------|
| Type             | -40 Deep Freezer |
| Cabinet Type     | Upright          |
| Climate Class    | N                |
| Cooling Type     | Direct cooling   |
| Defrost Mode     | Manual           |
| Refrigerant      | CFC-Free         |
| Noise Level (dB) | 45               |

### Performance

|                        |         |
|------------------------|---------|
| Cooling performance( ) | -40     |
| Temperature Range( )   | -20~-40 |

### Control

|            |                |
|------------|----------------|
| Controller | Microprocessor |
| Display    | LED            |

### Electrical Data

|                    |                       |
|--------------------|-----------------------|
| Power Supply(V/Hz) | 220-240/50/208-230/60 |
|--------------------|-----------------------|

|  |  |     |
|--|--|-----|
| <b>Power(W)</b>                          | 430                                    | 430 |
| <b>Electrical Current(A)</b>             | 2.9                                    | 2.9 |
| <b>Dimensions</b>                        |  |     |
| <b>Capacity(L/Cu.Ft)</b>                 | 490/17.3                               |     |
| <b>Net/Gross Weight(approx)</b>          | 137/145(kg)<br>374.8/440.9(lbs)        |     |
| <b>Interior Dimensions(W*D*H)</b>        | 685*650*1228(mm)<br>27.0*25.6*48.3(in) |     |
| <b>Exterior Dimensions(W*D*H)</b>        | 895*840*1860(mm)<br>35.2*33.1*73.2(in) |     |
| <b>Packing Dimensions(W*D*H)</b>         | 965*945*2040(mm)<br>38.0*37.2*80.3(in) |     |
| <b>Container load<br/>(20'/40'/40'H)</b> | 12/24/24                               |     |
| <b>Functions</b>                         |  |     |
| <b>High/Low Temperature</b>              | Y                                      |     |
| <b>Remote Alarm</b>                      | Y                                      |     |
| <b>Power Failure</b>                     | Y                                      |     |
| <b>Sensor Error</b>                      | Y                                      |     |
| <b>Low Battery</b>                       | Y                                      |     |
| <b>Accessories</b>                       |  |     |
| <b>Caster</b>                            | Y                                      |     |
| <b>Foot</b>                              | Y                                      |     |
| <b>Test Hole</b>                         | Y                                      |     |
| <b>Drawers</b>                           | 10                                     |     |
| <b>Temperature Recorder</b>              | Optional                               |     |
| <b>USB Interface</b>                     | N/A                                    |     |
| <b>Others</b>                            |  |     |
| <b>Certificate</b>                       | CE                                     |     |