



Key Features

- Reliable Performance
- Energy Saving
- Low Noise

More Detail:

http://www.haiermedical.com/products/preservation/40_deep_freezer/dw-40l626.html

Specifications

Technical Data

Type	-40 Deep Freezer
Cabinet Type	Upright
Climate Class	N
Cooling Type	Direct cooling
Defrost Mode	Manual
Refrigerant	CFC-Free
Noise Level (dB)	49

Performance

Cooling performance()	-40
Temperature Range()	-20~-40

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220-240/50
--------------------	------------

Power(W)	700
Electrical Current(A)	8
Dimensions	
Capacity(L/Cu.Ft)	626/22.1
Net/Gross Weight(approx)	298/314(kg) 657.0/692.2(lbs)
Interior Dimensions(W*D*H)	760*630*1310(mm) 29.9*24.8*51.6(in)
Exterior Dimensions(W*D*H)	1035*900*1980(mm) 40.7*35.4*78.0(in)
Packing Dimensions(W*D*H)	1072*954*2130(mm) 42.2*37.6*83.9(in)
Container load (20'/40'/40'H)	12/24/24
Functions	
High/Low Temperature	Y
Remote Alarm	Y
Hot Condenser	Y
Power Failure	Y
High/Low Voltage	Y
Sensor Error	Y
Low Battery	Y
High Ambient Temp	Y
Accessories	
Caster	Y
Foot	Y
Test Hole	Y
Shelves/Inner Doors	3/4
Temperature Recorder	Optional
USB Interface	N/A
Others	
Certificate	/