



Key Features

- Adjustable Temperature
- Five Alarm Functions
- Strong Refrigeration

More Detail:

http://www.haiermedical.com/products/preservation/40_deep_freezer/dw-40w100.html

Specifications

Technical Data

Type	-40 Deep Freezer
Cabinet Type	Chest
Climate Class	N
Cooling Type	Direct cooling
Defrost Mode	Manual
Refrigerant	CFC-Free

Performance

Cooling performance()	-40
Temperature Range()	-20~-40

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220-240/50
Power(W)	260

Electrical Current(A)	2.5
Dimensions	
Capacity(L/Cu.Ft)	100/3.5
Net/Gross Weight(approx)	43/46(kg) 94.8/101.4(lbs)
Interior Dimensions(W*D*H)	500*385*625(mm) 19.7*15.2*24.6(in)
Exterior Dimensions(W*D*H)	680*600*805(mm) 26.8*23.6*31.7(in)
Packing Dimensions(W*D*H)	735*645*855(mm) 28.9*25.4*33.7(in)
Container load (20'/40'/40'H)	54/108/162
Functions	
High/Low Temperature	Y
Remote Alarm	N/A
Hot Condenser	N/A
Power Failure	N/A
High/Low Voltage	N/A
Sensor Error	Y
Low Battery	N/A
High Ambient Temp	N/A
Accessories	
Caster	N/A
Foot	N/A
Test Hole	Y
Temperature Recorder	N/A
USB Interface	N/A
Others	
Certificate	CE