



Key Features

- Temperature Display And Alarm System
- Adjustable freezer compartment range: -20 ~ -40
- Adjustable Temperature

More Detail:

http://www.haiermedical.com/products/preservation/40_deep_freezer/dw-40w255.html

Specifications

Technical Data

Type	-40 Deep Freezer
Cabinet Type	Chest
Climate Class	N
Cooling Type	Direct cooling
Defrost Mode	Manual
Refrigerant	CFC-Free

Performance

Cooling performance()	-40
Temperature Range()	-20~-40

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220-240/50
Power(W)	410

Electrical Current(A)	2.5
Dimensions	
Capacity(L/Cu.Ft)	255/9.0
Net/Gross Weight(approx)	69/81(kg) 152.1/178.6(lbs)
Interior Dimensions(W*D*H)	1040*430*605(mm) 40.9*16.9*23.8(in)
Exterior Dimensions(W*D*H)	1243*633*838(mm) 48.9*24.9*33.0(in)
Packing Dimensions(W*D*H)	1310*750*880(mm) 51.6*29.5*34.6(in)
Container load (20'/40'/40'H)	24/54/81
Functions	
High/Low Temperature	Y
Remote Alarm	N/A
Hot Condenser	N/A
Power Failure	N/A
High/Low Voltage	N/A
Sensor Error	Y
Low Battery	N/A
High Ambient Temp	N/A
Accessories	
Caster	N/A
Foot	N/A
Test Hole	Y
Temperature Recorder	N/A
USB Interface	N/A
Others	
Certificate	CE