



Key Features

- Adjustable Temperature
- Temperature Display And Alarm System
- Five Alarm Functions

More Detail:

http://www.haiermedical.com/products/preservation/40_deep_freezer/dw-40w380.html

Specifications

Technical Data

Type	-40 Deep Freezer
Cabinet Type	Chest
Climate Class	N
Cooling Type	Direct cooling
Defrost Mode	Manual
Refrigerant	CFC-Free

Performance

Cooling performance()	-40
Temperature Range()	-20~-40

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220-240/50/208-230/60
Power(W)	550 550

Electrical Current(A)	3.4	3.4
Dimensions		
Capacity(L/Cu.Ft)	380/13.4	
Net/Gross Weight(approx)	82/89(kg) 180.8/196.2(lbs)	
Interior Dimensions(W*D*H)	1376*457*625(mm) 54.2*18.0*24.6(in)	
Exterior Dimensions(W*D*H)	1550*633*838(mm) 61.0*24.9*33.0(in)	
Packing Dimensions(W*D*H)	1630*755*880(mm) 64.2*29.7*34.6(in)	
Container load (20'/40'/40'H)	18/42/63	
Functions		
High/Low Temperature	Y	
Sensor Error	Y	
Accessories		
Caster	N/A	
Foot	N/A	
Test Hole	N/A	
Temperature Recorder	N/A	
USB Interface	N/A	
Others		
Certificate	CE	