



### Key Features

- Reliable Performance
- Energy Saving
- Slim Design

### More Detail:

[http://www.haiermedical.com/products/preservation/86\\_ult\\_freezer/dw-86l338.html](http://www.haiermedical.com/products/preservation/86_ult_freezer/dw-86l338.html)

## Specifications

### Technical Data

Type	-86 ULT Freeer
Cabinet Type	Upright
Climate Class	N
Cooling Type	Direct cooling
Defrost Mode	Manual
Refrigerant	CFC-Free
Noise Lever (dB)	50

### Performance

Cooling performance( )	-86
Temperature Range( )	-40~-86

### Control

Controller	Microprocessor
Display	LED

### Electrical Data

Power Supply(V/Hz)	220-240/50/208-230/60
--------------------	-----------------------

<b>Power(W)</b>	900	900
<b>Electrical Current(A)</b>	8	8
<b>Dimensions</b>		
<b>Capacity(L/Cu.Ft)</b>	338/11.9	
<b>Net/Gross Weight(approx)</b>	238/278(kg) 524.7/612.9(lbs)	
<b>Interior Dimensions(W*D*H)</b>	465*630*1165(mm) 18.3*24.8*45.9(in)	
<b>Exterior Dimensions(W*D*H)</b>	812*893*1846(mm) 32.0*35.2*72.7(in)	
<b>Packing Dimensions(W*D*H)</b>	875*970*1990(mm) 34.5*38.2*78.3(in)	
<b>Container load (20'/40'/40'H)</b>	12/24/24	
<b>Functions</b>		
<b>High/Low Temperature</b>	Y	
<b>Remote Alarm</b>	Y	
<b>Hot Condenser</b>	Y	
<b>Power Failure</b>	Y	
<b>High/Low Voltage</b>	Y	
<b>Sensor Error</b>	Y	
<b>Low Battery</b>	Y	
<b>High Ambient Temp</b>	Y	
<b>Door Ajar</b>	N/A	
<b>Power Consumption</b>	0.67Kw/h	
<b>Accessories</b>		
<b>Caster</b>	Y	
<b>Foot</b>	Y	
<b>Test Hole</b>	Y	
<b>Shelves/Inner Doors</b>	3/2	
<b>Temperature Recorder</b>	Optional	
<b>USB Interface</b>	N/A	
<b>Rs232/485 Interface</b>	Optional	
<b>CO2 Backup System</b>	Optional	

**LN2 Backup System**

Optional

**Racks And Boxes**

Optional

**Racks & Boxes(Optional)**

<b>Rack Model</b>	<b>Rack Description</b>	<b>Rack Dimensions (W*D*H)</b>	<b>Q'ty of Racks</b>	<b>Q'ty of Boxes</b>	<b>Q'ty of Samples 9*9 box / 10*10 box</b>
DCJ-54-A	the second and third layer:2 inch stainless steel rack 5*4	139*544*287(mm) 5.5*21.4*11.3(in)	6	120	9720/12000
DCJ-54-B	the second and third layer:2 inch stainless steel drawer rack 5*4	140*554*289(mm) 5.5*21.8*11.4(in)	6	120	9720/12000
DCJ-44-A	the first and forth layer:2 inch stainless steel rack4*4	139*544*234(mm) 5.5*21.4*9.2(in)	6	96	7776/9600
DCJ-44-B	the first and forth layer:2 inch stainless steel drawer rack4*4	140*554*236(mm) 5.5*21.8*9.3(in)	6	96	7776/9600

**Others****Certificate**

CE