



Key Features

- Reliable Performance
- Energy Saving
- Low Noise

More Detail:

http://www.haiermedical.com/products/preservation/86_ult_freezer/dw-86l486.html

Specifications

Technical Data

Type	-86 °C ULT Freezer
Cabinet Type	Upright
Climate Class	N
Cooling Type	Direct cooling
Defrost Mode	Manual
Refrigerant	CFC-Free
Noise Level (dB)	49

Performance

Cooling performance()	-86
Temperature Range()	-40~-86

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220-240/50/208-230/60/115/60
--------------------	------------------------------

Power(W)	1000	1000	1200
Electrical Current(A)	9	5.5	18
Dimensions			
Capacity(L/Cu.Ft)	486/17.1		
Net/Gross Weight(approx)	290/310(kg) 639.3/683.4(lbs)		
Interior Dimensions(W*D*H)	590*630*1310(mm) 23.2*24.8*51.6(in)		
Exterior Dimensions(W*D*H)	945*900*1980(mm) 37.2*35.4*78.0(in)		
Packing Dimensions(W*D*H)	990*985*2130(mm) 38.9*38.8*83.9(in)		
Container load (20'/40'/40'H)	12/24/24		
Functions			
High/Low Temperature	Y		
Remote Alarm	Y		
Hot Condenser	Y		
Power Failure	Y		
High/Low Voltage	N/A		
Sensor Error	Y		
Low Battery	Y		
High Ambient Temp	Y		
Door Ajar	Y		
Accessories			
Caster	Y		
Foot	Y		
Test Hole	Y		
Shelves/Inner Doors	3/4		
Temperature Recorder	Optional		
USB Interface	Y		
Rs232/485 Interface	Optional		
CO2 Backup System	Optional		
LN2 Backup System	Optional		

Racks And Boxes

Optional

Racks & Boxes(Optional)

Rack Model	Rack Description	Rack Dimensions (W*D*H)	Q'ty of Racks	Q'ty of Boxes	Q'ty of Samples 9*9 box / 10*10 box
DCJ-54-A	2 inch stainless steel rack 5*4	139*544*287(mm) 5.5*21.4*11.3(in)	16	320	25920/32000
DCJ-54-B	2 inch stainless steel drawer rack 5*4	140*554*289(mm) 5.5*21.8*11.4(in)	16	320	25920/32000

Others

Certificate

CE/UL