



Key Features

- Reliable Performance
- Energy Saving
- Low Noise

More Detail:

http://www.haiermedical.com/products/preservation/86_ult_freezer/dw-86l490.html

Specifications

Technical Data

Type	-86 ULT Freeer
Cabinet Type	Upright
Climate Class	N
Cooling Type	Direct cooling
Defrost Mode	Manual
Refrigerant	CFC-Free
Noise Lever (dB)	49

Performance

Cooling performance()	-86
Temperature Range()	-40~-86

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220-240/50/208-230/60
--------------------	-----------------------

Power(W)	1100	1100
Electrical Current(A)	6.5	6.5
Dimensions		
Capacity(L/Cu.Ft)	490/17.3	
Net/Gross Weight(approx)	295/335(kg) 650.4/738.5(lbs)	
Interior Dimensions(W*D*H)	590*630*1310(mm) 23.2*24.8*51.6(in)	
Exterior Dimensions(W*D*H)	860*900*1980(mm) 33.9*35.4*78.0(in)	
Packing Dimensions(W*D*H)	920*990*2130(mm) 36.2*38.9*83.9(in)	
Container load (20'/40'/40'H)	12/24/24	
Functions		
High/Low Temperature	Y	
Remote Alarm	Y	
Hot Condenser	Y	
Power Failure	Y	
High/Low Voltage	Y	
Sensor Error	Y	
Low Battery	Y	
High Ambient Temp	Y	
Door Ajar	N/A	
Accessories		
Caster	Y	
Foot	Y	
Test Hole	Y	
Shelves/Inner Doors	3/4	
Temperature Recorder	Optional	
USB Interface	N/A	
Rs232/485 Interface	Optional	
CO2 Backup System	Optional	
LN2 Backup System	Optional	

Racks And Boxes

Optional

Racks & Boxes(Optional)

Rack Model	Rack Description	Rack Dimensions (W*D*H)	Q'ty of Racks	Q'ty of Boxes	Q'ty of Samples 9*9 box / 10*10 box
DCJ-54-A	Upper chamber:2 inch stainless steel rack 5*4	139*544*287(mm) 5.5*21.4*11.3(in)	8	160	12960/16000
DCJ-54-B	Upper chamber:2 inch stainless steel drawer rack 5*4	140*554*289(mm) 5.5*21.8*11.4(in)	8	160	12960/16000
DCJ-44-A	Lower chamber:2 inch stainless steel rack4*4	139*544*234(mm) 5.4*21.4*9.2(in)	8	128	10368/12800
DCJ-44-B	the first and forth layer:2 inch stainless steel drawer rack4*4	140*554*236(mm) 5.5*21.8*9.3(in)	8	128	10368/12800

Others

Certificate

CE