



Key Features

- Touch Screen
- Low Noise
- Reliable Performance

More Detail:

http://www.haiermedical.com/products/preservation/86_ult_freezer/dw-86l729.html

Specifications

Technical Data

Type	-86 °C ULT Freezer
Cabinet Type	Upright
Climate Class	N
Cooling Type	Direct Cooling
Defrost Mode	Manual
Refrigerant	CFC-Free
Noise Level (dB)	49

Performance

Cooling performance()	-86
Temperature Range()	-40~-86

Control

Controller	Microprocessor
Display	LCD Touch Screen

Electrical Data

Power Supply(V/Hz)	220-240/50
--------------------	------------

Power(W)	1200
Electrical Current(A)	7
Dimensions	
Capacity(L/Cu.Ft)	729/25.7
Net/Gross Weight(approx)	310/332(kg) 683.4/731.9(lbs)
Interior Dimensions(W*D*H)	766*710*1310(mm) 27.0*28.0*51.6(in)
Exterior Dimensions(W*D*H)	1046*1000*1980(mm) 41.2*35.2*78.0(in)
Packing Dimensions(W*D*H)	1125*1110*2130(mm) 44.3*43.7*83.9(in)
Container load (20'/40'/40'H)	10/20/20
Functions	
High/Low Temperature	Y
Remote Alarm	Y
Hot Condenser	Y
Power Failure	Y
High/Low Voltage	Y
Sensor Error	Y
Low Battery	Y
High Ambient Temp	Y
Door Ajar	Y
Accessories	
Caster	Y
Foot	Y
Test Hole	Y
Shelves/Inner Doors	3/4
Temperature Recorder	Y
USB Interface	Y
Rs232/485 Interface	Optional
CO2 Backup System	Optional
LN2 Backup System	Optional

Racks And Boxes

Optional

Racks & Boxes(Optional)

Rack Model	Rack Description	Rack Dimensions (W*D*H)	Q'ty of Racks	Q'ty of Boxes	Q'ty of Samples 9*9 box / 10*10 box
DCJ-55-B	2 inch stainless steel drawer rack 5*5	140*688*289(mm) 5.5*27.1*11.4(in)	20	500	40500/50000
DCJ-55-A	2 inch stainless steel rack 5*5	138*680*287(mm) 5.4*26.8*11.3(in)	20	500	40500/50000

Others

Certificate

CE