



Key Features

- Touch Screen
- World Largest Capacity
- Reliable Performance

More Detail:

http://www.haiermedical.com/products/preservation/86_ult_freezer/dw-86l959.html

Specifications

Technical Data

Type	-86 ULT Freezer
Cabinet Type	Upright
Climate Class	N
Cooling Type	Direct Cooling
Defrost Mode	Manual
Refrigerant	CFC-Free
Noise Level (dB)	49

Performance

Cooling performance()	-86
Temperature Range()	-40~-86

Control

Controller	Microprocessor
Display	LCD Touch Screen

Electrical Data

Power Supply(V/Hz)	220-240/50
--------------------	------------

Power(W)	1300
Electrical Current(A)	9
Dimensions	
Capacity(L/Cu.Ft)	959/33.9
Net/Gross Weight(approx)	450/485(kg) 992.1/1069.2(lbs)
Interior Dimensions(W*D*H)	1016*710*1310(mm) 40.2*28.0*51.6(in)
Exterior Dimensions(W*D*H)	1296*998*1980(mm) 51.0*39.3*78.0(in)
Packing Dimensions(W*D*H)	1355*1100*2130(mm) 53.3*43.3*83.9(in)
Container load (20'/40'/40'H)	8/18/18
Functions	
High/Low Temperature	Y
Remote Alarm	Y
Hot Condenser	Y
Power Failure	Y
High/Low Voltage	Y
Sensor Error	Y
Low Battery	Y
High Ambient Temp	Y
Door Ajar	Y
Accessories	
Caster	Y
Foot	Y
Test Hole	Y
Shelves/Inner Doors	3/4
Temperature Recorder	Optional
USB Interface	Y
Rs232/485 Interface	Optional
CO2 Backup System	Optional
LN2 Backup System	Optional

Racks And Boxes

Optional

Racks & Boxes(Optional)

Rack Model	Rack Description	Rack Dimensions (W*D*H)	Q'ty of Racks	Q'ty of Boxes	Q'ty of Samples 9*9 box / 10*10 box
DCJ-115-A	2 inch stainless steel drawer rack 11*5	141*694*636(mm) 5.6*27.3*25.0(in)	14	770	62370/77000
DCJ-55-B	2 inch stainless steel drawer rack 5*5	140*688*289(mm) 5.5*27.1*11.4(in)	28	700	56700/70000
DCJ-55-A	2 inch stainless steel rack 5*5	138*680*287(mm) 5.4*26.8*11.3(in)	28	700	56700/70000

Others

Certificate

CE