



-86°C ULT Freezer DW-86W100



Key Features

- Reliable Performance
- Low Noise
- Easy To Use

More Detail:

http://www.haiermedical.com/products/preservation/86_ult_freezer/dw-86w100.html

Specifications

Technical Data

Type	-86 ULT Freezer
Cabinet Type	Chest
Climate Class	N
Cooling Type	Direct cooling
Defrost Mode	Manual
Refrigerant	CFC-Free
Noise Level (dB)	49

Performance

Cooling performance()	-86
Temperature Range()	-40~-86

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220-240/50/115/60	
Power(W)	650	650
Electrical Current(A)	7.5	12

Dimensions

--	--

Capacity(L/Cu.Ft)	100/3.5
Net/Gross Weight(approx)	138/160(kg) 304.2/352.7(lbs)
Interior Dimensions(W*D*H)	470*450*480(mm) 18.5*17.7*18.9(in)
Exterior Dimensions(W*D*H)	769*825*1120(mm) 30.3*32.5*44.1(in)
Packing Dimensions(W*D*H)	834*854*1235(mm) 32.9*33.9*48.6(in)
Container load (20'/40'/40'H)	12/24/48

Functions

High/Low Temperature	Y
Remote Alarm	Y
Hot Condenser	Y
Power Failure	Y
High/Low Voltage	Y
Sensor Error	Y
Low Battery	Y
High Ambient Temp	Y
Door Ajar	N/A

Accessories

Caster	Y
Foot	Y
Test Hole	Y
Shelves/Inner Doors	-/1
Temperature Recorder	Optional
USB Interface	N/A
Rs232/485 Interface	Optional
CO2 Backup System	Optional
LN2 Backup System	Optional
Racks And Boxes	Optional

Racks & Boxes(Optional)

Rack Model	Rack Description	Rack Dimensions (W*D*H)	Q'ty of Racks	Q'ty of Boxes	Q'ty of Samples 9*9 box / 10*10 box
-------------------	-------------------------	------------------------------------	----------------------	----------------------	--

DCJ-08-A	2 inch stainless steel rack 8*1	138*142*461(mm) 5.4*5.6*18.1(in)	9	72	5832/7200
----------	------------------------------------	-------------------------------------	---	----	-----------

Others

Certificate

CE/UL
