



Combined Refrigerator and Freezer HYCD-282A



Key Features

- Refrigerator & Freezer Combination
- Strong Refrigeration
- Five Alarm Functions

More Detail:

http://www.haiermedical.com/products/preservation/combined_refrigerator_and_freezer/hycd-282a.html

Specifications

Technical Data

Type	Combined Refrigerator and Freezer
Cabinet Type	Upright
Climate Class	N
Cooling Type	Refrigerator: Forced air cooling Freezer: Direct cooling
Defrost Mode	Manual/Auto
Refrigerant	CFC-Free
Noise Level (dB)	43

Performance

Temperature Range()	2~8/-20~-40
----------------------	-------------

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220-240/50/220/60
Power(W)	400 400

Electrical Current(A)	2.8	2.8
Dimensions		
Capacity(L/Cu.Ft)	Refrigerator: 185/6.5 Freezer: 97/3.4	
Net/Gross Weight(approx)	145/160(kg) 320.0/352.7(lbs)	
Interior Dimensions(W*D*H)	Refrigerator: 605*510*720(mm) 23.8*20.0*28.3(in) Freezer: 515*465*440(mm) 20.3*18.3*17.3(in)	
Exterior Dimensions(W*D*H)	736*660*1810(mm) 29.0*26.0*71.3(in)	
Packing Dimensions(W*D*H)	875*770*1920(mm) 34.4*30.3*75.6(in)	
Container load (20'/40'/40'H)	12/26/26	
Functions		
High/Low Temperature	Y	
Remote Alarm	Y	
Power Failure	Y	
Sensor Error	Y	
Low Battery	Y	
High Ambient Temp	Y	
Accessories		
Caster	Y	
Foot	Y	
Test Hole	Y	
Shelves/Inner Doors	3/-	
Drawers	2	
USB Interface	Y	
Others		
Certificate	/	